

# SRC Seating Jig

## Precision Seating for Screw-Retained Restorations

The SRC Seating Jig is a positioning index used to ensure accurate and efficient delivery of screw-retained crowns. Fabricated alongside the final restoration, it guides proper seating and orientation—reducing guesswork and chair time.

Used chairside, it helps:

- Ensure correct alignment with the implant connection
- Simplify delivery, especially in posterior or limited-visibility cases

**Made from visible light-cured material, it is provided by the lab with all SRC cases at no additional cost.**



## Why It Matters

### Faster Delivery

- Streamlines seating process
- Reduces chair time

### Improved Accuracy

- Ensures restoration seats as designed
- Minimizes risk of misalignment

### Simplified Workflow

- Especially beneficial in: Posterior regions, Limited access cases
- Provides a more controlled, predictable insertion